

myAIRVO 2 sets the standard for home humidification therapy.

With an integrated flow generator and (optional) oxygen delivery system, myAIRVO 2 Home based Humidification is an innovative system designed to enhance mucociliary clearance, optimise oxygen delivery and improve lung function.

As a high-performing solution for patients in the home environment and in long-term care, myAIRVO 2 is especially designed to support every breath you take.

With its integrated flow generator and innovative oxygen delivery system, the myAIRVO 2 Humidifier:

- Provides quality therapy for conditions such as COPD or bronchiectasis and anyone undergoing radiation therapy<sup>1</sup> where airway dryness or mucous retention is an issue
- Replicates hospital-grade treatment the high-performing myAIRVO 2 machine has been designed to exceed the ISO standard for respiratory humidification; ensuring maximum comfort and results for you in your home environment
- Delivers greater lifestyle freedom studies have shown that patients on long-term humidification therapy have significantly fewer exacerbation days, reduced exacerbation frequency and improved quality of life compared to usual care<sup>2</sup>

Call or email the CPAP Direct crew today to discuss how the myAIRVO 2 gives you the versatility and convenience to get on with your life.

1 – Andrew Macann, FRANZCR, et al., Phase 3 Trial of Domiciliary Humidification to Mitigate Acute Mucosal Toxicity During Radiation Therapy for Head-and-Neck Cancer: First Report of Trans Tasman Radiation Oncology Group (TROG) 07.03 RadioHUM Study

2 - Rea H, et al., The clinical utility of long-term humidification therapy in chronic airway disease, Respiratory Medicine (2010), doi:10.1016/j.rmed.2009.12.016

### myAIRVO 2 - Features



#### A complete system

 myAIRVO 2, with its wide range of accessories provides greater versatility, mobility and convenience



## Designed for simple set-up and use

- Simple controls
- Helpful on-screen animations assist with setup and troubleshooting



### **Versatile application**

- myAIRVO 2 can be used with Optiflow™ nasal interface as well as direct-connect tracheostomy and mask interfaces
- Wide flow range with automatic compensation for changes in interface
- Can add Oxygen from a cylinder or concentrator (when required)



# Adjustable temperature and flow settings

- Three temperature settings available (31°C, 34°C or 37°C) for added comfort and compliance
- Integrated flow generator delivers a wide flow range (2-60 L/min) for both paediatric and adult patients
- No wall air supply required

