



Thermoform 101.000

	102 Background White		100 White L		006 Metal Flake (± pms 873C)
	701 black opaque		700 Black		124/125 Silver (± pms 877C)
	200 Primrose Yellow (± pms 101C)		P200 Primrose Yellow (± pms Yellow C)		324 Pale Gold (± pms 873C)
	201 Lemon Yellow (± pms 108C)		P202 Chrome yellow (± pms 116C)		424 Rich Pale Gold (± pms 872C)
	202 Chrome Yellow (± pms 7548C)		P301 Orange (± pms 172C)		524 Rich gold (± pms 871C)
	301 Orange (± pms 1655C)		P302 Scarlet (± pms 185C)		824 Gold Flake (± pms 8643C)
	302 Scarlet (± pms Bright Red C)		P303 Red (± pms 185C)		850 Bronze Binder
	303 Red (± pms 199C)		P305 Carmine E (± pms 1935C)		880 Tixotropic Gel
	304 Ruby Red (± pms 187C)		P342 Pink (± pms 233C)		890 Overprint Varnish A.R.
	305 Carmine Red (± pms 188C)		P404 Bright Blue (± pms 286C)		
	342 Pink (± pms 675C)		P440 Violet (± pms 274C)		
	400 light Blue (± pms 2173C)		P500 green K (± pms 335C)		
	401 Cobalt Blue (± pms 300C)		215 Arc Chrome (± pms 804C)		
	404 Brilliant Blue (± pms 2945C)		310 Blaze		
	403 Ultramarine Blue (± pms 2736C)		312 Fire Orange (± pms 805C)		
	405 Dark Blue (± pms 294C)		317 Aurora Pink (± pms 806C)		
	440 Violet (± pms 275C)		511 Signal green (± pms 802C)		
	500 Leaf Green (± pms 3285C)		297 Tr. Yellow Tix (pms Process Yellow)		
	501 Light Green (± pms 354C)		397 Tr. Magenta Tix (pms Process Magenta)		
	503 Dark Green (± pms 3435C)		497 Tr. Cyan Tix (pms Process Cyan)		
	541 Green (± pms 348C)		797 Tr. Black Tix (pms Process Black)		