

CBRNe ORDNANCE & IED TRAINING KIT



Replica ordnance includes (from left-to-right): 120mm Mortar, 155mm Artillery Projectile, 122mm Rocket Warhead, 107mm Rocket, and 105mm Artillery Projectile



Inert improvised devices can be used as dispersal devices or as IEDs when combined with the included inert explosives and components


Next-Level Capabilities for CBRNe Response Training

- Chemical ordnance with integral containers for liquid simulants that can be implemented into "Leak, Seal, Pack" training scenarios
- Ordnance can be reconfigured into IEDs using the included inert explosives, detonators, & switches
- Improvised devices with a variety of switches can utilize liquid disperal or IED configurations, including time, victim, or command operation
- Each kit includes GB, HD, and VX Chemical Warfare Agent (CWA) simulants that provide an accurate detection response

MAC 7 and Inert Products work with personnel within various areas of conflict to bring forth the most up to date and complete training packages available today.



CWA Simulants



Unparalleled Realism

Familiarise detection teams with the visual and tactile properties of the most dangerous Chemical Warfare Agents and help responders recognise them easily during the course of their work.

This next generation of simulants not only gives identical results on the majority of detectors, but they also react accurately on detection paper.

Benefits Include:

- Identical look, feel, and behavior to the real chemicals
- Detection devices respond the same way as they would with real agent, such as Ion Mobility Spectrometers, Flame Spectrometers, or detection paper, eliminating the requirement to purchase expensive CBRN training detectors.
- Provides highly realistic training not normally affordable, but absolutely essential for today's professional.





EACH KIT CONTAINS:

A 250 ml bottle of GB, HD, and VX simulant A 20 ml dropper bottle for GB, HD, VX simulant A 10 ml demonstration Vial with GB, HD and VX

The 20 ml dropper bottles come in a separate user kit which can be used by the trainer to place out agent samples for their responders to discover and identify. They can be easily refilled from the 250 ml stock bottles.

The 10 ml clear glass demonstration vials enable the trainer to show responders what each liquid agent is like in terms of color, density and viscosity.

The CWA simulants are non-toxic when used by professional users, easy to use, and biodegradable.

MAC 7 and Inert Products work with personnel within various areas of conflict to bring forth the most up to date and complete training packages available today.