

The unique design of the Ultrasound Exam Step permits efficient and accurate exams of the upper and lower leg during venous insufficiency studies.

For the sonographer's comfort, the standing patient is placed in the optimal position for easy and ergonomically sound capture of the highest quality images. To switch legs the platform can be rotated 180 degrees using the nylon sliding feet.

- Designed for patients weighing up to 350 lbs / 158.7 kg
- Tip-N-Tow Easily portable with built-in rollers and grab bars
- Easy to keep clean, the platform is made from germ-resistant, non-porous polyethylene, and stainless steel

ULTRASOUND EXAM STEP - COMPACT
--------------------------------

Item #	Description	
24602	For standing leg venous insufficiency studies	
Dimens	1 <sup>st</sup> Step Height: 2 <sup>nd</sup> Step Height: Platform Width: Platform Depth:	14 in / 36 cm 24 in / 61 cm

*Clear Image Devices* has become the #1 brand in imaging accessories through innovative product development in cooperation with the medical imaging community. We believe that value begins with quality; quality products and customer service. We are advancing medical imaging every day by bridging the gap between technology and people.



Free Standard Shipping All Products - 48 States

Clear Image Devices <sup>LLC</sup> 5225 68th St., Suite 8 Caledonia, MI 49316 **Phone:** +1 734-474-6537 **Fax:** +1 734.661.0220



www.ClearImageDevices.com

copyright ©2016 Clear Image Devices LLC All rights reserved.