



[www.heyrex.com](http://www.heyrex.com)

**T.** +64 4 476 8849  
**F.** +64 4 476 2168  
**E.** [info@heyrex.com](mailto:info@heyrex.com)

**Designed in New Zealand**

Manufactured in China under Certified Supervision

**Patent Pending © 2006 Say International Ltd. All Rights Reserved**

Heyrex +

# TORUS™

*less stress, less mess*

**Product  
and  
Instruction  
Booklet**



# TORUS™

*less stress, less mess*



## To transport Torus™

6

**Tip:** After locking, remove any remaining water from the drinking well. On arrival, turn the Cap to the dog symbol.



Align the Cap arrow with the padlock symbol. This locks the water supply within the walls of **TORUS™** and prevents the flow of water into the drinking well.

(fig. 6)



5

## Torus™ Benefits

- + Large 2 Litre Water Capacity
- + Low, Stable Profile
- + Water Fill, Flow and Lock Options
- + Fresh Filtered Water
- + Made From Food Grade Materials
- + Reduces Water Spillage
- + Non-Slip Rubber Foot Pads
- + Ergonomic Shape Provides Easy Handling
- + Stored Water Protection
- + Convenient Water Portability



To release a flow of water into the drinking well, align the Cap arrow with the dog symbol.

(fig. 5)




# Why Torus™?

Dogs need plenty of fresh, clean water every day to remain healthy. As dog carers, it is our responsibility to ensure their needs are met.

Conventional water bowls often contain only a limited amount of water which can be reduced further by evaporation or spills. These problems decrease the water available and can lead to your dog being denied its daily drinking needs.

4



Replace the Cap by aligning the arrow with the tap symbol and inserting the Cap back into **TORUS™**.

(fig. 4)

Hold **TORUS™** horizontally and fill the reservoir with water through the opening.

3

**Tip:** Or you can fill **TORUS™** on the ground using a jug.



(fig. 3)

**TORUS™** (Latin for “hollow circular ring”) is the solution. Two litres of water stored inside the walls ensures an ample supply of water, protected from dust and contaminants. The low bowl profile also means no spills.

**TORUS™** is a “go-anywhere” bowl. The valve position, when locked, secures stored water inside the walls so you can take the filled **TORUS™** with you wherever you go.

**TORUS™ - less stress, less mess**



# How does Torus™ work?

The **TORUS™** design is comprised of a water storage area and a drinking well. The bulk of your dog's water is stored and protected inside the hollow walls of **TORUS™**. Each time your dog takes a drink, replacement water automatically flows from the storage area into the drinking well.

Activated carbon in the replaceable filter reduces flavours and trace elements commonly present in municipal water supplies.

**TORUS™** is exceptionally stable. Low profile geometry and soft rubber feet help prevent accidental spills by boisterous dogs or passers-by.

The water level in the drinking well is always high enough to allow easy drinking but low enough to minimise water spills and splashes.

2

**Tip:** Replacing the filter monthly is recommended.



Use a **TORUS™** replaceable filter to ensure your pet's water supply always remains fresh. Take a fresh filter and hold it for 10 seconds under a running tap to remove possible carbon fragments. Drop the replacement filter into place as shown.

(fig. 2)



# How to use your Torus™



Align the Cap arrow with the tap symbol and lift the Cap straight up to remove it.

(fig. 1)

1

1


Stored water is contained inside the **TORUS™** walls.

2

Fresh water is drawn from the walls through a carbon filter.

3

The filtered water flows out into the drinking well where it is kept at an optimum level.



# Colour Options

**TORUS™** is a very clever combination of simple, contemporary design elements and colour combinations that complement any home environment.

When buying a water bowl, always ask for **TORUS™**.

# What is Heyrex™?

Patented **Heyrex™** products are the foundation for our long term vision to provide premium pet care solutions and exclusive benefits.

With support from leading veterinary and University authorities around the globe, the name **Heyrex™** has become synonymous with personalised pet care for discerning pets.

