



Fully editable protocols - Perfect for OEM/ODM!

ZiXpress-32[®]


semi Automated Nucleic Acid Purification System

ZiXpress-32[®] is a compact and reliable automated extraction platform, based on the Magnetic Pillar technology that can extract up to 32 samples within 20 - 40 minutes with the optimized usage of the high-class *Zi-Beads*[®]. The flexible open system allows users to customize the protocols to meet precisely any application need.



- Best C/P value Extraction solution on the market
- Easy-To-Use User-Interface with 7" Touch Screen
- Advanced Magnetic Beads Extraction Technology
- Consistent and reproducible purified Results
- Cost-Effective and reliable

1 Cost Saving 	
	2 High Purity
3 Easy One Touch 	



Viral DNA/RNA



Bacterial



Whole Blood



Tissue

Features

Convenient & User-friendly	Ready-To-Go Solution with prefilled 96 well plate cartridges and consumables	✓
	Easy-To-Use User Interface multi-Language & intuitive UI with 7" Touchscreen - no external PC required	✓
	Ergonomic Instrument Design compact & light bench top instrument - fits in any laboratory	✓
Technology & Performance	Optimized Magnetic Pillar - Technology for reliable and reproducible extraction results	✓
	Perfect Match of Magnetic Particle Separation with Zinexts's customized Magnetic Beads <i>ZiBeads</i> ®	✓
	Customizable Open System Platform fully editable protocols for the use of your own Reagent Kits	✓
Low-Cost & Effectivity	Best Price and C/P Value get the most out of your money	✓



Applications




Infectious Disease



Genetic Testing



Oncology



Molecular Pathology



NIPT



Forensics



GMO & Food Safety



Quarantine

System Specifications

Extraction Method	Dimensions
Magnetic Beads Extraction Technology	40 cm W x 42 cm D x 44cm H

Extraction Steps	Weight
Sample binding, washing and elution	32 kg

Throughput	Operating Temperature
1 to 32 samples per run	0-30°C

Display	Operating Relative Humidity
7" LCD Touchpanel	30-80 % (non-condensing)

Heat Block Temperature	Electrical Requirements
Up to 200°C	110-240V 50/60 Hz

Prefilled Extraction Kits (each 64 tests)

ZiXpress-32® Viral DNA & RNA Extraction Kit	ZP02201
ZiXpress-32® Bacteria Genomic DNA Extraction Kit	ZP02202
ZiXpress-32® Whole Blood Genomic DNA Extraction Kit	ZP02203
ZiXpress-32® Tissue Genomic DNA Extraction Kit	ZP02204



<https://twitter.com/zinexts>
<https://linkedin.com/company/zinexts>
<https://facebook.com/zinexts>