

QUICK NITRATE & NITRITE TEST SAMPLING INSTRUCTIONS

Quick Nitrate & Nitrite is a Quick Series test package designed to provide rapid test results.

PLEASE READ ALL INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR WATER SAMPLE

✓ This test kit represents one sampling location. If you need to test multiple locations, you will need to purchase multiple test kits.

Samples must be shipped the same day they are collected via guaranteed next day delivery

DO NOT collect and/or ship your water samples on Friday, Saturday or the day before a National Holiday.

 \checkmark To ensure valid bacteria results, overnight shipping arrangements should be made prior to beginning your water sample collection.

Samples should be collected after 1:00 p.m. Eastern Time.

- 1. Freeze the blue gel ice pack for at least 8 hours prior to sample collection. Return to the cooler just prior to sampling and/or return shipping.
- 2. Collect samples on Monday through Thursday, as late in the day as possible before you ship the package for overnight delivery. Samples *MUST* arrive at our laboratory within 24 hours of collection.
- 3. Remove aerators, strainers, mixers or purifiers from your tap. Disinfect spigot using alcohol or flame (if metal). *If you suspect methane in your water do not use flame.*
- 4. Run the cold water for approximately 3-5 minutes, then adjust the flow to a thin stream to avoid splashing or overfilling.
- 5. Fill the 500 ml bottle to the shoulder, cap securely and return it to the cooler.
- 6. Complete the enclosed lab form, **making sure to indicate the collection date and time you sampled**, type of water (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Retain a copy of the lab form for your records and return the completed form to the shipping container.
- 7. Mail your sample to the testing location below. The sample should be sent via guaranteed overnight delivery.

Don't forget to complete and return the lab form with the sample shipment! Save money on shipping by purchasing at reduced rate. Call NTL at 800-458-3330 for pricing.

> Testing Location: National Testing Laboratories, Ltd. 556 South Mansfield Street Ypsilanti, MI 48197

8. PLEASE NOTE: Testing typically takes approximately 2-3 business days from receipt of samples. Results will be emailed upon completion.

NOTE: The lab reserves the right to reject water samples that appear to be excessively contaminated as they may cause potential damage to laboratory equipment.

If you have any questions regarding this Quick Nitrate & Nitrite Test or any other services we provide, please contact a Technical Service Representative at National Testing Laboratories at 1-800-458-3330 for assistance.

Thank you for choosing National Testing Laboratories, Ltd., for your water testing services!

National Testing Laboratories, Ltd. 1-800-H2O-TEST / 1-800-458-3330 www.ntllabs.com / www.watercheck.com ntlsales@ntllabs.com

