

Handpiece

- Ergonomic design, autoclavable sleeve

Sterilized, Disposable Fiber Tips for Every Procedure

- 200 μ m (21mm), 300 μ m (9mm), 400 μ m (9mm) diameters

Display Stand

- Mount on Handpiece rack for easy and convenient viewing

Wireless Foot Pedal

- Easy-to-use, convenient placement

Temporary Pain Relief*

- Relieve maxillofacial or intraoral discomfort
- Reduce inflammation
- Three on-device presets for biostimulation, inflammation, and muscle/joint relief
- Relief is immediate

Paired with an iPod Touch

- iPod touch included
- 21 precise procedure presets
- Set-up guide and tutorials
- Training and procedure videos include 8CE credits, a \$795 value
- Unlimited upgrade capability
- Contact CAO customer service via phone, email, and Facetime® on the iPod touch
- Clinical support from certified diode laser professionals
- Utilize iPod touch technology in your office for clinical photos, videos, professional applications, etc.

Base

- Light weight, portable between operatories
- iPod touch charger w/USB port
- Bluetooth Connection
- Mountable on dental chair
- Powerful performance with 3 Watts of power for all soft-tissue procedures

Warranty

- Two years parts and labor*
See warranty card for details.

**Extended warranty available.
Please contact us.*

Precise White™ Laser Assisted In-Office Teeth Whitening

- Whiten up to 5 shades in one 30 minute procedure
- 30% hydrogen peroxide gel
- Pairs with Sheer White! for a complete whitening system
- Can be focused on a single tooth or the full arch
- Precise Dam™ included in Precise White kit
- Can be used for intracoronal bleaching

