


# PREMIUM Palm Protector with cylinder roll




The effective moisture & pressure reduction solution.



- Product provides immediate moisture control as well as cushioning support inside the hand
- Cylinder roll will prevent the hand from closing completely thereby improving flexion contractures inside the hand.
- Product will keep the skin at a regulated temperature which will eliminate the likelihood of skin breakdown
- Product will eliminate odors generally caused by excess humidity or contractures
- Dri-release fibers found inside the product effectively attract and expel moisture 4 times faster than regular fabric
- The 3-dimensional structure inside the product allows for air circulation inside the product
- Fabric Content: 90% Polyester 10% Cotton
- Washing temperature 160 – 167°F (71 – 75 °C)  
\* Drying Temperature 220°F

LATEX FREE

Recommended for treatment of contracted hand



<input type="checkbox"/>	Right Hand
<input type="checkbox"/>	Left Hand

## Size Chart

- |                          |  |
|--------------------------|--|
| <input type="checkbox"/> | <b>SMALL</b> (3-4" Knuckle to knuckle / .75" Diameter / 4" Length)   |
| <input type="checkbox"/> | <b>MEDIUM</b> (4-5" Knuckle to knuckle / 1" Diameter / 4" Length)    |
| <input type="checkbox"/> | <b>LARGE</b> (5-6" Knuckle to knuckle / 1.5" Diameter / 4.5" Length) |

**TexMedico**  
www.texmedico.com

T: 1-506-855-9959 • F: 1-506-382-0079  
Toll Free: 1.877.448.4432 (CANADA/USA) Dieppe, NB, Canada